EN - English		1
FR - Français		14
ES - Español		28
IT - Italiano		42
DE - Deutsch		56
PT - Português		70
NL - Nederlands		84
NO - Norsk		89
DA - Dansk		94
SV - Svenska		99
FI - Suomi		104
HU - Magyar		109
PL - Polski		114

Revision 2
24/09/13

Instruction Manual

Rocking Platform Shakers

3-D Rotator Wavers

Incubating Rocking Platform Shakers

Incubating 3-D Rotator Wavers

European Catalogue Numbers:

Rocker Platform Shaker

Euro Plug: 444-0756

UK Plug: 444-0757

Swiss Plug: 444-0758

3-D Rotator Wavers

Euro Plug: 444-0759

UK Plug: 444-0760

Swiss Plug: 444-0761

Incubating Rocking Platform Shaker

Euro Plug: 444-0762

UK Plug: 444-0763

Swiss Plug: 444-0764

Incubating 3-D Rotator Waver

Euro Plug: 444-0765

UK Plug: 444-0766

Swiss Plug: 444-0767

TABLE OF CONTENTS

Package Contents	1
Warranty	1
Installation	1
Maintenance & Servicing	2
Intended Use	2
Equipment Disposal	2
Environmental Conditions	2
Safety Instructions	3
Warning	3
Standards & Regulations	3
Control Panel	4
Incubating Control Panel	5
Rocker Specifications	6
Waver Specifications	6
Incubating Rocker Specifications	7
Incubating Waver Specifications	7
Rocker/Waver Operating Instructions	8-9
Incubating Rocker/Waver Operating Instructions	10-11
Technical Service	12
Troubleshooting	12
Accessories	13
EC Declaration of Conformity	119-121

PACKAGE CONTENTS

- Rocking Platform Shaker or 3-D Rotator Waver or Incubating Rocking Platform Shaker or Incubating 3-D Rotator Waver
- 234cm detachable power cord
- Non-skid rubber mat (not included with incubating units)
- Instruction Manual

WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

INSTALLATION

Upon receiving the VWR Rocker/Waver, check to ensure that no damage has occurred in shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the Rocker/Waver on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of six (6) inches from vertical surfaces. Always place the unit on a sturdy work surface.

The Rocker/Waver is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The 230V unit plugs into a 230 volt, 50/60 Hz source.

Legal Manufacturer:

VWR International bvba • Researchpark Haasrode 2020 • Geldenaakesbaan 464 • B-3001 Leuven • + 32 16 385011 • <http://be.vwr.com>

MAINTENANCE & SERVICING

The Rocker/Waver is built for long, trouble-free, dependable service. No lubrication or other technical user maintenance is required. It needs no user maintenance beyond keeping the surfaces clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly. Do not use a cleaning agent or solvent on the front panel which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

INTENDED USE

These Rockers and Wavers are intended for general laboratory use.

EQUIPMENT DISPOSAL

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

ENVIRONMENTAL CONDITIONS - OPEN AIR ROCKERS AND WEVERS

Operating Conditions: Indoor use only.

Temperature*:	For use in CO ₂ environments, incubators or cold rooms ranging from -10 to 60°C
Humidity:	80% relative humidity, non-condensing
Altitude:	0 to 2000 M above sea level

Non-Operating Storage:

Temperature*:	-20 to 65°C
Humidity:	maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

* **Avoid cold starts:** Unit is not designed to start after being in a cold room environment. Bring unit into cold room from a room temperature environment, operate and remove unit from cold room as soon as operation is complete.

ENVIRONMENTAL CONDITIONS - INCUBATING ROCKERS AND WEVERS

Operating Conditions: Indoor use only.

Temperature:	5 to 40°C
Humidity:	80% relative humidity, non-condensing
Altitude:	0 to 2000 M above sea level

Non-Operating Storage:

Temperature:	-20 to 65°C
Humidity:	maximum 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

SAFETY INSTRUCTIONS

Please read the entire instruction manual before operating the Rocker/Waver.

WARNING! DO NOT use the Rocker/Waver in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the platform tray or lid.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly. Biohazard spills should be cleaned up using approved laboratory procedures. Solvent spills are a fire hazard. Stop the unit immediately, and **DO NOT** operate until clean up is complete and vapors have dissipated (motor brushes arc and may ignite flammable vapors).

DO NOT immerse the unit for cleaning.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground - Protective Conductor Terminal

Alternating Current

STANDARDS & REGULATIONS

VWR International hereby declares under its sole responsibility that the products comply with the requirements of the following directives and associated standards:

Associated EU guidelines:

EMC directive 2004/108/EC

LVD directive 2006/95/EC

Safety standards:

IEC 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use. Part I: General Requirements.

IEC 61010-2-010* Part II: Particular requirements for laboratory equipment for the heating of materials.

IEC 61010-2-051 Part II: Particular requirements for laboratory equipment for mixing and stirring.

UL Std. No. 61010-1

CSA/CAN C22.2 No. 0-M91

CSA/CAN C22.2 No. 61010-1-04

EMC standards:


EN61000-3-2 EN61000-4-3

EN61000-3-3 EN61000-4-4

EN61000-4-5 EN61000-4-6

EN61000-4-2 EN61000-4-11

* Applies to Incubating Rockers/Wavers only


ROCKER/WAVER CONTROL PANEL

The front panel of the Rocker/Waver contains all the controls and displays needed to operate the unit.

- A. Standby button/standby indicator light:** The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to activate the speed/tilt and time functions. The standby indicator light will shut off and the speed/tilt display and time display will illuminate. Press the standby button again and the unit will once again be in standby mode.
- B. Speed/tilt display:** Displays the speed and tilt of the unit. **C.** Up/down arrows for set-point control. On/off button starts/stops the rocking/waving function.
- D.** The speed/tilt indicator light will be illuminated when the unit is rocking/waving.

- E. Function button:** Press to choose the function you are setting: speed or tilt.
- F. Speed indicator light:** Illuminates when the speed is displayed.
- G. Tilt indicator light:** Illuminates when the tilt angel is displayed.
- H. Time display:** Displays accumulated time (continuous mode) or how much time is remaining (timed mode). **I.** Up/down arrows for set-point control. On/off button starts/stops the time function. The display range is from 0 to 9999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999.

INCUBATING ROCKER/WAVER CONTROL PANEL

The front panel of the Incubating Rocker/Waver contains all the controls and displays needed to operate the unit.

- A. Standby button/standby indicator light:** The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to activate the temperature, speed/tilt and time functions. The standby indicator light will shut off and the temperature, speed/tilt and time displays will illuminate. Press the standby button again and the unit will once again be in standby mode.
- B. Temperature display:** Displays the actual/set-point temperatures in conjunction with the actual/set-point indicator lights. **C.** Up/down arrows for set-point control. On/off button starts/stops the heat function. **D.** The heat indicator light will be illuminated when the unit is heating.
- E. Actual indicator light:** Illuminates when the temperature displayed is the actual temperature of the air in the chamber.

- F. Set-point indicator light:** Illuminates when the set-point temperature is displayed.
- G. Caution hot indicator light:** Illuminates when the air temperature of the chamber is above 40°C (104°F).
- H. Speed/tilt display:** Displays the speed and tilt of the unit. **I.** Up/down arrows for set-point control. On/off button starts/stops the rocking/waving function. **J.** The speed/tilt indicator light will be illuminated when the unit is rocking/waving.
- K. Function button:** Press to choose the function you are setting: speed or tilt.
- L. Speed indicator light:** Illuminates when the speed is displayed.
- M. Tilt indicator light:** Illuminates when the tilt angle is displayed.
- N. Time display:** Displays accumulated time (continuous mode) or how much time is remaining (timed mode). **O.** Up/down arrows for set-point control. On/off button starts/stops the time function. The display range is from 0 to 9999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999.

ROCKER SPECIFICATIONS

Overall dimensions (L x W x H):	43,2 x 27,9 x 12,7cm
Platform dimensions (L x W):	32,4 x 25,4cm
Electrical (50/60 Hz):	230 volts AC, 5 amps, 25 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Speed range:	1 to 50rpm*
Speed accuracy:	±1rpm
Tilt angle:	0 to 15°*
Weight capacity:	4,5kg**
Timer:	digital, 1 second to 9999 minutes (increased in 1 second increments)
Controls:	see page 4
Ship weight:	7kg
European Catalogue Numbers:	Euro plug: 444-0756
	UK plug: 444-0757
	Swiss plug: 444-0758

* Maximum speed/tilt angle may vary with heavy or unbalanced loads.

** Centered on tray.

WAVER SPECIFICATIONS

Overall dimensions (L x W x H):	43,2 x 27,9 x 15,2cm
Platform dimensions (L x W):	29,9 x 22,2cm
Electrical (50/60 Hz):	230 volts AC, 5 amps, 25 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Speed range:	1 to 30rpm*
Speed accuracy:	±1rpm
Tilt angle:	0 to 20°*
Weight capacity:	2,3kg**
Timer:	digital, 1 second to 9999 minutes (increased in 1 second increments)
Controls:	see page 4
Ship weight:	7.3kg
European Catalogue Numbers:	Euro plug: 444-0759
	UK plug: 444-0760
	Swiss plug: 444-0761

* Maximum speed/tilt angle may vary with heavy or unbalanced loads.

** Centered on tray.

INCUBATING ROCKER SPECIFICATIONS

Overall dimensions (L x W x H):	43,2 x 27,9 x 26,7cm
Platform dimensions (L x W):	25,4 x 19,1cm
Electrical (50/60 Hz):	230 volts AC, 5 amps, 450 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Temperature range:	ambient +5° to 65°C
Temperature uniformity:	+/-0,5°C at 37°C
Speed range:	1 to 50rpm*
Speed accuracy:	±1rpm
Tilt angle:	0 to 15°**
Weight capacity:	4,5kg**
Timer:	digital, 1 second to 9999 minutes (increased in 1 second increments)
Controls:	see page 5
Ship weight:	10kg
European Catalogue Numbers:	Euro plug: 444-0762 UK plug: 444-0763 Swiss plug: 444-0764

* Maximum speed/tilt angle may vary with heavy or unbalanced loads.

** Centered on tray.

INCUBATING WAVER SPECIFICATIONS

Overall dimensions (L x W x H):	43,2 x 27,9 x 26,7cm
Platform dimensions (L x W):	23,5 x 18,4cm
Electrical (50/60 Hz):	230 volts AC, 5 amps, 450 watts
Fuses:	5mm x 20mm, 5 amp quick acting
Temperature range:	ambient +5° to 65°C
Temperature uniformity:	+/-0,5°C at 37°C
Speed range:	1 to 30rpm*
Speed accuracy:	±1rpm
Tilt angle:	0 to 20°*
Weight capacity:	2,3kg**
Timer:	digital, 1 second to 9999 minutes (increased in 1 second increments)
Controls:	see page 5
Ship weight:	10kg
European Catalogue Numbers:	Euro plug: 444-0765 UK plug: 444-0766 Swiss plug: 444-0767

* Maximum speed/tilt angle may vary with heavy or unbalanced loads.

** Centered on tray.

ROCKER/WAVER OPERATING INSTRUCTIONS

The VWR Rocker/Wavers are used for gentle mixing of laboratory samples. These units have been designed for the speed/tilt and time functions to work independently of one another. The speed/tilt can be re-set without re-setting the timer and the timer can be stopped and started without interrupting the rocking/waving functions.

1. Getting ready:

- a. Plug the power cord into a properly grounded outlet. The standby indicator light will illuminate, verifying power to the unit. Manually move the platform tray to a near horizontal "home" position before powering unit.
- b. Press the standby button to move the unit from standby mode. The standby indicator light will shut off and the speed/tilt and time displays will illuminate and display the previously used settings.

2. Setting tilt: Electronic tilt enables the user to adjust the angle of tray movement while the unit is rocking/waving or while the unit is stopped.

- a. Press the function button below the speed/tilt display until the tilt indicator light illuminates. You are now ready to set the tilt angle.
- b. Press the up/down arrows below the speed/tilt display until you reach the desired angle. When you release the arrow button, the display will blink off and then on indicating the new set angle has been accepted. The unit will complete one rotation at the previously set angle, then smoothly change to the newly programmed angle. The tilt indicator light will flash until the transition to the new tilt angle is complete.

3. Setting speed:

- a. Press the function button below the speed/tilt display until the speed indicator light illuminates. You are now ready to set the speed.
- b. Press the up/down arrows below the speed/tilt display until you reach the desired speed. When you release the arrow button, the display will blink off and then on indicating the new set speed has been accepted.

- c. Press the on/off button to start the rocking/waving function. The indicator light below the speed/tilt display will flash until the transition to the set speed is complete, then the light will stay illuminated indicating the rocking/waving function is in use.
- d. To stop rocking/waving function, press the on/off button below the speed/tilt display. The unit will complete one full rotation then stop in the horizontal "home" position. The speed indicator light will flash until the rotation is complete. When the cycle is complete, the unit will automatically move to standby mode and the standby indicator light will illuminate.

OPERATING TIPS

When the unit is running at slow speeds or high angles, making large changes to the tilt angle or speed may take several minutes to complete. The quickest way to make large changes to the speed or tilt angle is to stop the unit, change the speed or tilt angle, then restart by pressing the on/off button below the speed/tilt display.

4. Setting time to zero (0:00) and continuous mode: Accumulated time.

- a. Press and hold the on/off button below the time display. After three (3) seconds, the display will indicate the previous set time.
- b. Simultaneously press both the up and down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00) minutes. Alternately, you can use the up/down arrows to get to zero (0:00).
- c. Press the on/off button below the time display. The display will indicate the accumulated time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will **NOT** interrupt the rocking/waving function. Press the on/off button below the speed/tilt display to interrupt the rocking/waving function.
- d. To re-set, press and hold the on/off button below the time display. After three (3) seconds, the display will indicate the previous set time, which was zero (0:00).

ROCKER/WAVER OPERATING INSTRUCTIONS (CONT'D)

5. Setting timed mode: Programmed time.

- a. Press the up/down arrows below the time display until you reach the desired time.
- b. Start this function by pressing the on/off button below the time display. The unit will run for the selected time, the up/down arrows will become inactive while the timer is running. The unit will stop rocking/waving when time display reaches zero (0:00). Four (4) audible beeps will indicate the time down function is complete. The time display will default back to the set time. To repeat for the same time, simply depress the on/off button again.
- c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the time display. The time display will flash off and on to indicate the time function is on "hold". **IMPORTANT:** This will **NOT** interrupt the rocking/waving function. Press the on/off button below the speed/tilt display to interrupt the rocking/waving function. Restart the timer by pressing the on/off button below the time display. Unit will continue counting down to zero (0:00). When the display reaches zero (0:00), you will hear the four (4) audible beeps that indicate the time down function is complete and the rocking/waving function will cease.

6. Turning unit off:

- a. To turn the unit off, press the standby button. The speed/tilt and time displays will be blank, the standby indicator light will illuminate. The Rocker/Waver should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord from the unit or un-plug from the wall outlet.

OPERATING TIPS

Built-in memory maintains the last used speed/tilt and time settings during a power interruption.

A built-in program will shut power off to the motor if the platform tray is prevented from rocking/waving, or the unit is overloaded beyond it's recommended weight capacity.

INCUBATING ROCKER/WAVER OPERATING INSTRUCTIONS

The VWR Incubating Rocker/Wavers are used to prepare samples for testing. These units have been designed for the temperature, speed/tilt and time functions to work independently of one another. The temperature and speed/tilt can be re-set without re-setting the timer, and the timer can be stopped and started without interrupting the heating or rocking/waving functions.

1. Getting ready:

- a. Plug the power cord into a properly grounded outlet. The standby indicator light will illuminate, verifying power to the unit. Manually move the platform tray to a near horizontal "home" position before powering unit.
- b. Press the standby button to move the unit from standby mode. The standby indicator light will shut off and the temperature, speed/tilt and time displays will illuminate and display the previously used settings.

2. Set temperature:

- a. Press the up/down arrows below the temperature display until you reach the desired temperature. When you release the arrow button, the display will blink off and then on indicating the new set temperature has been accepted.
- b. Press the on/off button to start the heating function. The indicator light below the temperature display will illuminate to indicate the heating function is in use. You will hear five (5) audible beeps that indicate set temperature has been reached.
- c. Set-point temperature adjustments can be made without interrupting heating using the up/down arrows below the temperature display. After the change has been made and you release the button, the display will blink off and then on indicating the new set temperature has been accepted.
- d. To stop heating, press the on/off button below the temperature display, the heat indicator light will turn off.

Caution Hot Indicator:

The caution hot indicator light warns that the temperature of the air in the chamber is above 40°C. The light will illuminate and remain on when the temperature of the air in the chamber reaches approximately 40°C. When the heat is turned off, the caution hot indicator light will stay illuminated until the temperature of the air in the chamber is less than 40°C.

3. **Setting tilt:** Electronic tilt enables the user to adjust the angle of tray movement while the unit is rocking/waving or while the unit is stopped.

- a. Press the function button below the speed/tilt display until the tilt indicator light illuminates. You are now ready to set the tilt angle.
- b. Press the up/down arrows below the speed/tilt display until you reach the desired angle. When you release the arrow button, the display will blink off and then on indicating the new set angle has been accepted. The unit will complete one rotation at the previously set angle, then smoothly change to the newly programmed angle. The tilt indicator light will flash until the transition to the new tilt angle is complete.

4. Setting speed:

- a. Press the function button below the speed/tilt display until the speed indicator light illuminates. You are now ready to set the speed.
- b. Press the up/down arrows below the speed/tilt display until you reach the desired speed. When you release the arrow button, the display will blink off and then on indicating the new set speed has been accepted.
- c. Press the on/off button to start the rocking/waving function. The indicator light below the speed/tilt display will flash until the transition to the set speed is complete, then the light will stay illuminated indicating the rocking/waving function is in use.

INCUBATING ROCKER/WAVER OPERATING INSTRUCTIONS (CONT'D)

- d. To stop rocking/waving function, press the on/off button below the speed/tilt display. The unit will complete one full rotation then stop in the horizontal "home" position. The speed indicator light will flash until the rotation is complete. When the cycle is complete, the unit will automatically move to standby mode and the standby indicator light will illuminate.

OPERATING TIPS

When the unit is running at slow speeds or high angles, making large changes to the tilt angle or speed may take several minutes to complete. The quickest way to make large changes to the speed or tilt angle is to stop the unit, change the speed or tilt angle, then restart by pressing the on/off button below the speed/tilt display.

5. **Setting time to zero (0:00) and continuous mode:** Accumulated time.
 - a. Press and hold the on/off button below the time display. After three (3) seconds, the display will indicate the previous set time.
 - b. Simultaneously press both the up and down arrows, the display will indicate zero (0:00). The unit time is now set to zero (0:00) minutes. Alternately, you can use the up/down arrows to get to zero (0:00).
 - c. Press the on/off button below the time display. The display will indicate the accumulated time. The up/down arrows will become inactive. To stop timer, press the on/off button again. **IMPORTANT:** This will **NOT** interrupt the rocking/waving function. Press the on/off button below the speed/tilt display to interrupt the rocking/waving function.
 - d. To re-set, press and hold the on/off button below the time display. After three (3) seconds, the display will indicate the previous set time, which was zero (0:00).
6. **Setting timed mode:** Programmed time.
 - a. Press the up/down arrows below the time display until you reach the desired time.

- b. Start this function by pressing the on/off button below the time display. The unit will run for the selected time, the up/down arrows will become inactive while the timer is running. The unit will stop rocking/waving when time display reaches zero (0:00). Four (4) audible beeps will indicate the time down function is complete. The time display will default back to the set time. To repeat for the same time, simply depress the on/off button again.
- c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the time display. The time display will flash off and on to indicate the time function is on "hold". **IMPORTANT:** This will **NOT** interrupt the rocking/waving function. Press the on/off button below the speed/tilt display to interrupt the rocking/waving function. Restart the timer by pressing the on/off button below the time display, unit will continue counting down to zero (0:00). When the display reaches zero (0:00), you will hear the four (4) audible beeps that indicate the time down function is complete and the rocking/waving function will cease.

7. Turning unit off:

- a. To turn unit off, press the standby button. The temperature, speed/tilt and time displays will be blank, the standby indicator light will illuminate. The Incubating Rocker/Waver should be kept in standby mode when not in use. To completely cut off power to the unit, disconnect the power cord from the unit or un-plug from the wall outlet.

OPERATING TIPS

Built-in memory maintains the last used temperature, speed/tilt and time settings during a power interruption.

A built-in program will shut power off to the motor if the platform tray is prevented from rocking/waving, or the unit is overloaded beyond its recommended weight capacity.

TECHNICAL SERVICE

For more information or technical assistance contact your local VWR representative or visit www.vwr.com.

Web resources: Visit VWR's website at www.vwr.com for:

- Complete technical service contact information
- Access to VWR's Online Catalogue, and information about accessories and related products
- Additional product information and special offers

TROUBLESHOOTING

During operation, any rattling or ticking sounds may indicate a loose screw on the platform tray, tray attachment or accessory. All accessories should be sufficiently tightened in place before starting the unit.

Error	Cause of Error	How to Fix
E1	RTD open or temperature over 100°C for an Incubating Unit	This error cannot be fixed by the end user. Unplug the unit and please contact your VWR representative for repair.
E2	RTD shorted or temperature below 0°C for an Incubating Unit	This error cannot be fixed by the end user. Unplug the unit and please contact your VWR representative for repair.
E3	Mechanical obstruction Ceased motor Drive belt broken	Press the standby button to clear this error. If the E3 error persists, along with grinding, knocking or rubbing noises, un-plug the unit and contact your VWR representative for repairs. In the event of an errant E3 error, such as someone inadvertently touching the platform tray while the unit is running, the unit will automatically go to the horizontal "home" position and restart. In the event of an obstruction left under the tray, which would cause the unit to continually restart, the unit will attempt to auto home and restart four (4) times and then stop running, displaying an E3 error. This error can be cleared by the user by pressing the standby button.
E4	Maximum load exceeded	Press the standby button to clear this error. Be sure the load is within the maximum load specification before restarting the unit. If problem persists, unplug the unit and please contact your VWR representative for repair.


ACCESSORIES

DIMPLED MAT

Designed to hold centrifuge tubes, vials, culture tubes, and micro-tubes securely in place. Mat can be easily removed for cleaning and transporting tubes from bench to tray.

DESCRIPTION	DIMENSIONS	PART NUMBER
Dimpled Mat for Rocking Platform Shaker	32,4 x 25,4 cm	444-0768
Dimpled Mat for 3-D Rotator Waver	29,9 x 22,2 cm	444-0770
Dimpled Mat for Incubating Rocking Platform Shaker	25,4 x 19,1 cm	444-0769
Dimpled Mat for Incubating 3-D Rotator Waver	23,5 x 18,4 cm	444-0771

STACKING TRAY *(not for use with Incubating Rockers or Wavers)*

Easily attaches to the lower platform tray to add a second tier for higher capacity applications. Second tier tray mounts 8,9 cm above the lower tray.

DESCRIPTION	PART NUMBER
Stacking Tray for Rocking Platform Shaker	444-0772
Stacking Tray for 3-D Rotator Waver	444-0773

EN - English		1
FR - Français		14
ES - Español		28
IT - Italiano		42
DE - Deutsch		56
PT - Português		70
NL - Nederlands		84
NO - Norsk		89
DA - Dansk		94
SV - Svenska		99
FI - Suomi		104
HU - Magyar		109
PL - Polski		114

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hereby declare that the products, mentioned in Schedule 1 comply with the essential requirements of the following Directives of the European Parliament and council:

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EN61010-1:2001	Relating to safety requirements for electrical equipment for measurement, control and laboratory use—part 1: general requirements— second edition, 2001, including amendment 1 (1992) and amendment 2 (1995)
IEC 61010-2-010:2005	Relating to safety requirements for electrical equipment for measurement, control and laboratory use—Part 2-010: particular requirements for laboratory equipment for heating of materials

1/3

VWR International bvba/sprl Researchpark Haasrode 2020, Geldenaaksaan 464, B-3001 Leuven, Belgium
 Tel.: 016 385 011 • Fax: 016 385 385 • E-mail: info@be.vwr.com

Handwritten signature and date:
 16 April 2013

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IEC 61010-2-051:2005	Relating to safety requirements for electrical equipment for measurement, control and laboratory use— part 2-051: particular requirements for laboratory equipment for mixing and stirring
IEC 61000-3-2, IEC 61000-3-3, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5, IEC 61000-4-6, IEC 61000-4-1	Relating to electromagnetic compatibility

Year of CE marking: **2006**

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

Place and date of issue: **22.04.2013**

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Yves Van Damme
European Regulatory Affairs Manager

Crisina Silva
Managing Director

EC Declaration of Conformity

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Schedule 1
Annexe 1
Anhang 1

Article Number	Article description
444-0756	ROCKER- ROCKING PLATFORM SHAKER 230V EU ROCKER- ROCKING PLATFORM SHAKER 230V EU ROCKER- ROCKING PLATFORM SHAKER 230V EU
444-0757	ROCKER- ROCKING PLATFORM SHAKER 230V UK ROCKER- ROCKING PLATFORM SHAKER 230V UK ROCKER- ROCKING PLATFORM SHAKER 230V UK
444-0758	ROCKER- ROCKING PLATFORM SHAKER 230V CH ROCKER- ROCKING PLATFORM SHAKER 230V CH ROCKER- ROCKING PLATFORM SHAKER 230V CH
444-0759	WAVER- 3-D ROTATOR WAVER 230V EU WAVER- 3-D ROTATOR WAVER 230V EU WAVER- 3-D ROTATOR WAVER 230V EU
444-0760	WAVER- 3-D ROTATOR WAVER 230V UK WAVER- 3-D ROTATOR WAVER 230V UK WAVER- 3-D ROTATOR WAVER 230V UK
444-0761	WAVER- 3-D ROTATOR WAVER 230V CH WAVER- 3-D ROTATOR WAVER 230V CH WAVER- 3-D ROTATOR WAVER 230V CH
444-0762	ROCKER- INCUB PLATFORM SHAKER 230V EU ROCKER- INCUB PLATFORM SHAKER 230V EU ROCKER- INCUB PLATFORM SHAKER 230V EU
444-0763	ROCKER- INCUB PLATFORM SHAKER 230V UK ROCKER- INCUB PLATFORM SHAKER 230V UK ROCKER- INCUB PLATFORM SHAKER 230V UK
444-0764	ROCKER- INCUB PLATFORM SHAKER 230V CH ROCKER- INCUB PLATFORM SHAKER 230V CH ROCKER- INCUB PLATFORM SHAKER 230V CH
444-0765	WAVER INCUB- 3-D ROTATOR WAVER 230V EU WAVER INCUB- 3-D ROTATOR WAVER 230V EU WAVER INCUB- 3-D ROTATOR WAVER 230V EU
444-0766	WAVER INCUB- 3-D ROTATOR WAVER 230V UK WAVER INCUB- 3-D ROTATOR WAVER 230V UK WAVER INCUB- 3-D ROTATOR WAVER 230V UK
444-0767	WAVER INCUB- 3-D ROTATOR WAVER 230V CH WAVER INCUB- 3-D ROTATOR WAVER 230V CH WAVER INCUB- 3-D ROTATOR WAVER 230V CH

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Ireland / Northern Ireland

VWR International Ltd / VWR International
(Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: 01 88 22 222
Fax: 01 88 22 333
E-mail: sales@ie.vwr.com

Italy

VWR International PBI S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: 02-3320311/02-487791
Fax: 800 152999/02-40090010
E-mail: info@it.vwr.com
info@internationalpbi.it

The Netherlands

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: 020 4808 400
Fax: 020 4808 480
E-mail: info@nl.vwr.com

Norway

VWR International AS
Haavard Martinsens vei 30
0978 Oslo
Tel.: 02290
Fax: 815 00 940
E-mail: info@no.vwr.com

Poland

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: 058 32 38 200 do 204
Fax: 058 32 38 205
E-mail: labart@pl.vwr.com

Portugal

VWR International - Material de Laboratório,
Lda
Edifício Neopark
Av. Tomás Ribeiro, 43- 3 D
2790-221 Camaxide
Tel.: 21 3600 770
Fax: 21 3600 798/9
E-mail: info@pt.vwr.com

Spain

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Linars Park
08450 - Linars del Vallès
Barcelona
Tel.: 902 222 897
Fax: 902 430 657
E-mail: info@es.vwr.com

Sweden

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: 08 621 34 00
Fax: 08 621 34 66
E-mail: kundservice@se.vwr.com

Switzerland

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: 044 745 13 13
Fax: 044 745 13 10
E-mail: info@ch.vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: 0800 22 33 44
Fax: 01455 55 85 86
E-mail: uksales@uk.vwr.com

China

VWR (Shanghai) Co., Ltd
2nd Floor, Building 4,
Lane 998, Halei Rd,
Zhangjiang Hi-tech Park
Shanghai, 201203
China
Tel.: +86-21 589 868 88
Fax: +86-21 585 588 01
E-mail: info_china@vwr.com

India

VWR Lab Products Private Limited
135/12, Brigade Towers, 2nd Floor
Front wing, Brigade Road,
Bengaluru, India - 560 025
Tel.: +91-80-41117125/26 (Bengaluru)
Tel.: +91-2522-647911/922 (Mumbai)
Fax: +91-80-41117120
E-mail: vwr_india@vwr.com

Singapore

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
Fax: +65 6264 3780
E-mail: sales@sg.vwr.com

Legal Manufacturer:

VWR International bvba • Researchpark Haasrode 2020 • Geldenaaksebaan 464 • B-3001 Leuven • + 32 16 385011 • <http://be.vwr.com>
715105-00 (Rev 2)