

Datasheet

HEXIM1 monoclonal antibody (M01), clone 4G4

Catalog Number: H00010614-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant HEXIM1.

Clone Name: 4G4

Immunogen: HEXIM1 (AAH06460, 1 a.a. ~ 359 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MAEPFLSEYQHQPQTSNCTGAAAVQEELNPERPPGA
EERVPEEDSRWQSRAFPQLGGRPGPEGESLESQPP
PLQTQACPESSCLREGEKGQNGDDSSAGGDFPPPAE
VEPTPEAELLAQPCHDSEASKLGAPAAGGEEEWGQQ
QRQLGKKKHRRRPSKKKRHWKPYKLTWEEKKKFDE
KQSLRASRIRAEMFAKGQPVAPYNTTQFLMDDHDQEE
PDLKTGLYSKRAAAKSDDTSDDDFMEEGGEEDGGSD
GMGGDGSEFLQRDFSETYERYHTESLQNMSKQELIKE
YLELEKCLSRMEDENNRLRLESKRLGGDDARVRELEL
ELDRLRAENLQLLTENELHRQQERAPLSKFGD

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10614

Gene Symbol: HEXIM1

Gene Alias: CLP1, EDG1, FLJ13562, HIS1, MAQ1

Gene Summary: Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth muscle cells. This gene has no introns. [provided by RefSeq]