

Datasheet

SPAG5 monoclonal antibody (M01), clone 5F9

RefSeq]

Catalog Number: H00010615-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant SPAG5.

Clone Name: 5F9

Immunogen: SPAG5 (AAH00322, 1082 a.a. ~ 1191 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

AETETKVLQEALAGQLDSNCQPMATNWIQEKVWLSQ
EVDKLRVMFLEMKNEKEKLMIKFQSHRNILEENLRRSD
KELEKLDDIVQHIYKTLLSIPEVVRGCKELQGLLEF

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG3 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10615

Gene Symbol: SPAG5

Gene Alias: DEEPEST, MAP126, hMAP126

Gene Summary: This gene encodes a protein
associated with the mitotic spindle apparatus. The
encoded protein may be involved in the functional and
dynamic regulation of mitotic spindles. [provided by