



Product Specification Sheet

VWR FILTER TIPS RACKED

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Type	Capacity (µl)	Length(mm)	Inner Diameter (mm)	Packed	Pk
732-0542	UltraFine™, extended	0.1 - 10	38.1	0.25	10 racks, 96 each	960
732-0515	UltraFine™, graduated	0.1 - 10	31.4	0.51	12 racks, 96 each	1.152
732-0518	Bevel Point™, graduated	0.5 - 10	44.8	0.51	10 racks, 96 each	960
732-0527	Bevel Point™	1 - 40	49.6	0.51	10 racks, 96 each	960
732-0568	Bevel Point™	1 - 50	49.6	0.55	10 racks, 96 each	960
732-0522	Bevel Point™	1 - 100	49.6	0.55	10 racks, 96 each	960
732-0529	Bevel Point™	1 - 160	49.6	0.51	10 racks, 96 each	960
732-0559	Bevel Point™, wide orifice	1 - 200	49.6	1.45	10 racks, 96 each	960
732-0540	UltraFine™	1 - 200	52	0.51	10 racks, 96 each	960
732-0536	Bevel Point™	1 - 300	52	0.51	10 racks, 96 each	960
732-0533	UltraFine™	100 - 1000	88.9	0.76	6 racks , 96 each	576

Material:

- Polypropylene made of 100% virgin resin
- The resins used to produce these products have been tested for heavy metals contamination using the U.S.P. (231) method and found to have results lower than 1ppm.



Product Specification Sheet

Other Characteristics:

- Made in the highest quality moulds, these tips can be used whenever the work demands a higher standard.
- A wider seal assures that tips will seal easily and consistently without leaking, on old as well as new pipettes
- Autoclavable

Clean Claims:

- Certified testing is performed lot by lot.
- Free of detectable Nuclease contaminants (RNase, DNase) as confirmed by a gel assay
- Free from Human DNA to a detection level of 1 pico-gram of DNA
- Certificates are on file at the production facility.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **January 2014**