

Material 84893.290
Material description Silica gel 60 - 200 µm

CAS Number 7631-86-9
Molecular mass 60.08

Characteristics	Specifications
Bulk density (tapped)	460 --> 550 g/l
pH (5 %)	6.5 --> 7.5
Loss on drying (160°C; 15 min)	2.0 --> 4.0 %
Particle size (< 63 µm)	Max. 5.0 %
Particle size (> 200 µm)	Max. 5.0 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Anja Vanhalle, Head of Laboratory - Haasrode
VWR International bvba; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.