




SLC5A3 Antibody

CATALOG NUMBER: 48-167



Immunohistochemistry staining of SLC5A3
in spinal cord tissue using SLC5A3
Antibody.

Specifications

SPECIES REACTIVITY:	Gorilla, Human
TESTED APPLICATIONS:	IHC
APPLICATIONS:	SLC5A3 antibody can be used in ELISA, Western Blot starting at 1:500 - 1:1000, immunohistochemistry starting at 5 ug/mL, and immunofluorescence starting at 10 ug/mL.
USER NOTE:	Optimal dilutions for each application to be determined by the researcher.
SPECIFICITY:	BLAST analysis of the peptide immunogen showed no homology with other human proteins.
IMMUNOGEN:	SLC5A3 antibody was raised against a peptide located near the internal domain of SLC5A3 (Human).
HOST SPECIES:	Rabbit

Properties

PURIFICATION:	Immunoaffinity Chromatography
PHYSICAL STATE:	Liquid
BUFFER:	PBS, 0.1% sodium azide.
STORAGE CONDITIONS:	SLC5A3 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
CLONALITY:	Polyclonal
CONJUGATE:	Unconjugated

Additional Info

ALTERNATE NAMES:	SLC5A3, Na ⁺ /myo-inositol cotransporter, SMIT1, SMIT2, SMIT, BCW2 protein
ACCESSION NO.:	P53794
PROTEIN GI NO.:	3334498
OFFICIAL SYMBOL:	SLC5A3
GENE ID:	6526

Background

